



Bioluminescent whole-cell biosensors for point-of-need applications: new analytical devices exploiting 3D printing technology


Luca Cevenini, Maria Maddalena Calabretta, Elisa Michelini, Aldo Roda

Laboratory of Analytical and Bioanalytical Chemistry

Dip. Chimica G.Ciamician, University of Bologna

Istituto Nazionale Biostrutture e Biosistemi, INBB.

Bioluminescent whole-cell bioreporters



Peculiarities

- **Easy to engineer** to respond to different analytes
- Analyte **bioavailability**
- **Biological activity** of a sample
- **Low cost**, self replicating
- **Multiplexing** capability


Main issues

- Genetic stability of engineered cells (**low robustness**)
- Long assay and **response time**
- Limited **dynamic range** for quantitative analysis


Bioluminescence as detection principle

Bioluminescent reactions:

- high quantum yield emission → **high detectability**
- high specific signal (low background) → **high S/N ratio**
- no excitation light source is required → **simple instrumentation**
- different BL systems & spectral tunable emission → **multiplexing**



e.g. Firefly Luciferase catalyzed BL reaction



λ_{\max} : 560 nm, quantum yield: ~ 0.44



Analytical issue: development of portable BL biosensors

From lab...



Miniaturization

... to on-site analysis



Robustness

Sensitivity

Response time

- Bioreporters integration with a detection system
- Cell movement/control on chip

3D Printing Technology



Bioprinting: living tissue, organs



Guns



Bionic medical devices

Is 3D printing the future of everything?



Food



Houses



Space station spare parts

3D-Printing Analytical Devices

Desktop 3D printer




Makerbot Replicator 2X

- Fused polymer deposition
- Dual extrusion
- Resolution: 100µm

Filament (1.75mm):

- ABS (1kg, 25€)
- PLA
- Rubber



Makerware software

Molds for PDMS casting



~0,10 €
15min

Microfluidic chips



~0,05 €
10min

Multiwell cartridges




~0,15 €
20min




Development of magnetic/BL whole-cell bioreporters

Magnetotactic bacteria (MTB) are microorganisms ($0.5 \times 5 \mu\text{m}$) that can orient along geomagnetic field lines thanks to a **magnetosome chain**.

The magnetosome consists of **magnetite crystal** (Fe_3O_4 ; 30-50nm) surrounded by **lipid membrane**




Magnetic control of BL-MTB






Cells can be easily moved/trapped inside a microfluidic chip by applying external magnetic field

3D-printed Magnetox chip prototype

A **3D printed** microfluidic chip is created using the dual-extrusion mode: **black** and **transparent ABS** are fused to create a chip prototype.



BL-MTB can be moved to detection areas thanks to neodymium-iron-boron (NdFeB) magnets



[Roda A. et. al, *Lab Chip*, 2013]

Integration with a portable detection system

ATIK 11000

Sensor Type: CCD - Kodak 11002

Sensor size: 37.25mm x 25.70mm



Resolution: 4008 pixels x 2672 pixels

Pixel Size: 9 μM x 9 μM

Cooling: Peltier with $\Delta T = -38^\circ\text{C}$



Magnetox biosensor prototype



Bioluminescence detection with Smartphones

Advantages


- Large diffusion & Connectivity
- Applications for data acquisition/handling

Limitations

- Field of view, macro focusing
- No cooling system
- Short exposure time (few sec)



Samsung Galaxy S2



POOR SENSITIVITY


How to increase sensitivity: new BL reporters

Cloning and mutagenesis
of new BL genes



Synthetic luciferases

- Different emission
- Increased stability
- pH insensitive emission




Rational design



Luciferase	λ_{\max} (nm)	$T_{1/2}$ (h, 37°C)
PpyWT	557	0,3
LucGREEN	548	> 5
LucRED	615	> 5



Generation of BL whole-cell bioreporters for smartphone

Constitutive expression vector



Transient transfections


Hek293T cell line



Toxicity test with smartphone

Method

- Cells seeding (20 μ l) in the multiwell cartridge
- Incubation (30min) with increasing concentration of ZnSO₄ (5 μ l)
- Substrate addition: 10 μ l D-luciferin
- Signal acquisition: 4s
- Image elaboration : ImageJ
- Toxicity curve





Perspectives

Analytical applications of smartphone-based whole-cell biosensors




Heavy metal detection



Water analysis



Anti-doping test



Food contamination



Endocrine disruptors
Xenobiotics detection